SAA-CHECK-1 ad us. vet.

(Für/for EasyReader® / EasyReader Plus®)

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1.1. TEST-KIT COMPONENTS

1 test-kit **SAA-CHECK-1** contains:

1. 10 or 20 test cassettes, coated with monoclonal antibodies
2. 10 or 20 dropper bottles **A** with 2.0 ml buffer diluent each
3. 5 or 10 dropper bottles **A** (without buffer diluent)
4. 20 or 40 pipette tips (10 μ l)
5. 5 or 10 pipette tips (200 μ l)
6. 1 instructions for use

1.2. STABILITY AND STORAGE

Store at
15–25°C

Expiry date
– see label

For veterinary use only

Lot number

In vitro diagnosticum

Do not use test-kit components from different kits, lot numbers or beyond stated expiry date

Follow instructions for use carefully

1.3. LIABILITY

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from the use of this product.

2. INTRODUCTION

As response to tissue injury and inflammation, the body produces, among others, so-called acute-phase proteins (APP) in the liver. The Serum Amyloid A (SAA) is an important APP in many species, with low serum concentrations during normal homeostasis and a quick answer during and after inflammation.

The physiological SAA concentration in the cat is below 10 mg/l, in the horse below 30 mg/l. A drastic increase in the SAA level up to over 1000 times occurs within few hours (4–6 hours, peak 36 to 48 hours later) due to inflammations, infections, traumata or wound healing disorders (horse). As soon as the inflammation subsides through a successful therapy, the SAA concentration decreases (half-life 24–48 hours).

The veterinary SAA measurement in blood allows not only the determination of clinical or subclinical inflammations, but also a reliable monitoring of course and therapy. Another advantage of the SAA determination is its independence to antibiotic, glucocorticoids or stress.

3. TEST PRINCIPLE

SAA-CHECK-1 is a rapid test for the quantitative determination of Serum Amyloid A in whole blood, plasma or serum of the cat and horse.

By the combination of mobile conjugate-bound antibodies and membrane-bound antibodies, the SAA concentration in the sample is determined with a high degree of sensitivity.

Depending on SAA concentration in the sample, a test and a control line of different intensity will appear. They can be evaluated quantitatively **only with EasyReader®**.

4. INFORMATION ON THE SAMPLE MATERIAL



- The sample material should be collected under the standard laboratory conditions (aseptic, non-hemolytic, lipaemic or icteric).
 - **SAA-CHECK-1** can be performed with anticoagulated (EDTA, citrate or heparin) whole blood, plasma or serum of the dog. **Native blood (whole blood without anticoagulant) should not be used!**
 - **ATTENTION:** Partially filled and/or insufficient mixed EDTA, Citrate or Heparin tubes could create interference effects (overdose of anticoagulant) resulting in un-specific test results.
 - **Whole blood (with anticoagulant) must be tested within 4 hours.**
 - **ATTENTION:** At **15–25°C**, plasma or serum can be stored up to 24 hours. If the test is to be run after 24 hours, the sample should be stored at minimum **–20°C**. **Do not store in the fridge!** Avoid repeated freezing and thawing.
 - Mix the sample material well before use!
 - Keep in mind that the sample material, as well as all used test-kit components, should have reached room temperature (15–25°C) at the time of application.
-

5. EASYREADER® PREPARATION

- Before starting with sample preparation/test procedure, please read the Instructions for use of the Reader! The Reader can be used in **Plus Mode (PM)** and in **Easy Mode (EM)**: see chapter 12).
- For switching from EM to PM, press the blue **Plus Mode** button in the upper right of the display.
- For switching from PM to EM, first press the **Main** button, then the **Logout** button, then in the bottom left the blue **Immediate** button.
- **Important: Switch the Reader on and allow it to warm up for at least 30 minutes! The Reader is preset so that the automatic swith-off time is deactivated to guarantee a constant incubation temperature in the Reader.** If you wish a defined switch-off time, you can change it as described in chapter 11.10 of the Reader's instructions for use.
- Have the **blank device** ready! In an emergency, a test cassette freshly removed from the foil pouch can also be used.

6. SAMPLE PREPARATION

1. Label one dropper bottle containing buffer diluent with animal's ID.
2. Unscrew the dropper bottle.
- 3a. **Plasma / serum:** Add 10 μl of sample with help of the laboratory pipette into the dropper bottle (fig.1).
- 3b. **Whole blood (with anticoagulant):** Add 10 μl of sample with help of the laboratory pipette into the dropper bottle (fig.1). Repeat this step with a new pipette tip.
4. Replace the screw cap on the dropper bottle.
5. Mix the sample-buffer mixture (SBM) **strongly** by inverting the dropper bottle upside-down several times for 10 seconds (fig.2).


7. TEST PROCEDURE

7.1. EM

1. Choose the desired test type (here SAA).
2. Insert the blank cassette as required and remove it after successful blank measurement.
3. Enter patient's data into the Reader.
4. Choose **Countdown** or **Immediate**.
 - Countdown:** The internal timer of the EasyReader starts automatically. The test cassette is read after 5 minutes incubation time. **Countdown** is recommended, because the EasyReader has a constant incubation temperature inside.
 - Immediate:** You have to stop the incubation time yourself with a stop watch. As soon as the time is ready, place the test cassette into the test cassette holder and start the immediate measurement (**Start**).
5. Remove the test cassette from the foil pouch shortly before use.
6. Label the test cassette with the ID of the animal.
7. Break the tip of the dropper bottle containing the SBM (fig.3)
8. Hold the dropper bottle **vertically** and slowly add 4 drops (150 μl) of SBM into the sample window of the test cassette (fig.4). Avoid bubbles and touching the test cassette!
9. **Countdown:** As soon as the pink-purple lateral flow is visible in the result window (fig.5), place the test cassette, with sample window on the right, immediately into the test cassette holder (fig.6) and place this into the Reader. Incubation time starts running automatically.
 - Immediate:** As soon as the pink-purple lateral flow is visible in the result window (fig.5), start your stop watch with the indicated incubation time (5 minutes). After

the incubation time, place the test cassette, with sample window on the right, immediately into the test cassette holder (fig.6) and place this into the Reader. Attention! The EasyReader immediately starts the measurement!

- The test result (in mg/l) can be read after 5 minutes directly on the Reader and is printed automatically after finishing the incubation time.

7.2. PM (is automatically in Countdown by default)

- Choose the desired test type (here SAA).
- Insert the blank cassette as required and remove it after successful blank measurement.
- Remove the test cassette from the foil pouch shortly before use.
- Label the test cassette with the ID of the animal.
- Break the tip of the dropper bottle containing the SBM (fig.3)
- Hold the dropper bottle **vertically** and slowly add 4 drops (150 μ l) of SBM into the sample window of the test cassette (fig.4). Avoid bubbles and touching the test cassette!
- As soon as the pink-purple lateral flow is visible in the result window (fig.5), place the test cassette, with sample window on the right, immediately into the test cassette holder (fig.6), place this into the Reader and press **Start**. Incubation time starts running automatically.
- Enter patient's data into the Reader.
- The test result (in mg/l) can be read after 5 minutes directly on the Reader and is printed automatically after finishing the incubation time.

7.3. SERIAL TESTING OF SEVERAL SAMPLES

- Several tests can be scheduled at the same time. However, it is recommended that you test the same parameter for each test run.
- Each test cassette is started one minute after the other. The whole reaction time as well as the delay between the single test cassettes must be followed exactly (stop watch!).

- The test result must be measured exactly at the end of the incubation time to get a reliable test result.

7.3.1 EM

- For the testing of several samples, it is recommended to use an external timer (stop watch). **Immediate** must be activated on the Reader (see instructions for the Reader).

7.3.2 PM

- For serial testing in PM, first you have to activate and confirm **Immediate** in the menu Settings » QC Options. Then, in the measurement menu, tick the box **Immediate**.

8. READING OF THE TEST RESULT

The test result (in mg/l) can be read after 5 minutes directly on the Reader and is printed automatically after finishing the incubation time.

The quantitative measuring range is 5–75 mg/l. The SAA concentrations reported outside this range are summarised in ranges (see table 8.1).

SAA concentration	Reader results
0–5 mg/l	< 5 mg/l
5–75 mg/l	Quantitative results
75–150 mg/l	75–150 mg/l
> 150 mg/l	> 150 mg/l

Table 8.1

9. FURTHER DILUTION WITH RESULT ABOVE 75 MG/L


If the SAA concentration of the sample is above 75 mg/l, a second measurement with a diluted sample is recommended to return to the quantitative range. The result of the diluted sample must be calculated using the respective dilution factor (see table 9.2).

Reader result	dilution stage	addition of buffer diluent into empty drop-per bottle	SBM	final dilution	Test result Reader and calculation
< 5 mg/l	not necessary				
5–75 mg/l	not necessary				
75–150 mg/l	1:3 (corresponds to 1+2)	400 µl	200 µl	1:600	x mg/l × 3
> 150 mg/l	1:5 (corresponds to 1+4)	800 µl	200 µl	1:1000	x mg/l × 5

Table 9.2

Sample preparation of dilution

1. Label one **empty** dropper bottle with animal's ID.
2. Unscrew the **empty** dropper bottle.
3. For a 1:3 dilution, take a new, buffer diluent-filled dropper bottle. Open it and add 400 μl of buffer diluent into the empty dropper bottle (fig.7). For a 1:5 dilution, add 800 μl of buffer diluent into the empty dropper bottle.
4. Take the same pipette tip and remove 200 μl from the original sample-buffer mixture (SBM). Add these 200 μl into the dropper bottle already filled with 400 μl (1:3 dilution) or 800 μl (1:5 dilution) buffer diluent (fig.8).
5. Replace the screw cap on the dropper bottle.
6. Mix the SBM **strongly** by inverting upside-down the bottle several times for 10 seconds (fig.9).
7. Perform the test as described in chapter 7 "test procedure". Multiply the specified result with the respective dilution factor (see table 9.2).

10. INFORMATION FOR THE INTERPRETATION

Cat

- SAA level below 10 mg/l → diseases/disorders are unlikely (1, 2).
- Elevated SAA levels:
 - Injury: mean value 62 mg/l
 - Renal failure: mean value 31 mg/l
 - Infectious diseases: mean value 48 mg/l
 - Tumours: mean value 17 mg/l

Horse

- SAA level below 30 mg/l → systemic inflammation is unlikely
- SAA level within 37 and 1,609 mg/l (mean value: 343 mg/l) → local inflammation is probable.
- SAA level within 688 and 4,000 mg/l (mean value: 1,583 mg/l) → systemic inflammation is developing and requires attention (3).

11. FEATURES OF PERFORMANCE

a) Accuracy

A study has been performed using pre-assayed SAA canine samples by serial dilution in a SAA depleted serum and covering a SAA range of 0–150 mg/l.

Optical densities expressed as a function of SAA concentrations are described by following curve: $Y = 0.0002x^3 - 0.0832x^2 + 10.929x - 13.591$ ($r^2 = 1$)

b) Sensitivity

The **SAA-CHECK-1** test is detecting SAA concentration of 5 mg/l.

c) Precision

37 cat serum samples and 30 horse serum samples were assayed using the ROCHE Cobas analyser and **SAA-CHECK-1** quantitative test. The obtained results show a correlation coefficient of 94.8% (cat samples) and 87.4% (horse samples).

d) Hook effect

Samples having a SAA concentration up to 600 mg/l constantly showed "> 150 mg/l" results.

e) Intra assay reproducibility

Intralot precision was evaluated running 25 replicates of 3 pre-assayed serum samples containing 13.4 mg/l, 26 mg/l and 58 mg/l of SAA. The obtained CV's (coefficient of variation) are 12.7%, 8.1% and 6.7%, respectively.

12. PRECAUTIONS

- Label sample material and associated dropper bottle and test cassette to ensure a precise assignment.
- Use a new pipette tip and a new dropper bottle for each sample.
- The buffer diluent contains low concentrations of toxic sodium azide as a preservative, therefore avoid skin/eye contact and/or ingestion.
- The sample material must be seen as potentially infectious and disposed of accordingly, together with the used test-kit components.
- Wear protective clothing such as laboratory coats and disposable gloves while assaying samples. Wash the hands after finishing the test procedure.
- Do not eat, drink or smoke in the area where samples and kit reagents are handled.
- Avoid any contact between hands and eyes or nose during sample collection and testing.

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
* anticoagulated whole blood (EDTA, heparin, citrate)

* anti koaguleret helblod (EDTA, heparin, citrat)