


CopperBody[®]

Provides a barrier to contamination & patient interference to support a rapid recovery



Ergonomic Design
with stretch Fit



No Cone = Full freedom
of movement



30 UVA + UVB
Protection



Easy-On
Step-In Zipper with
Coat Protector



4-Legs
Design



Castration
Design



Embedded copper
protects & promotes
healthy skin



Washable 40°C



Breathable



Optional
reinforced
opening for
males

LUNA & MILO[®]

LOVING PET CARE

CopperBody®

Because they deserve it!



Copper Body Classic

Our original model suitable for males and females with optional opening for male dogs



Copper Body Male Dog Castration

Ideal for castration of male dogs as the design gives extra protection around the castration site



Copper Body 4 Legs

Ideal for full body protection as well as difficult to protect leg wounds or skin disorders


SIZE CHART

Size	Weight*	Back Length	Chest Circ.	Classic	Castration	4-Legs
00	0 - 1,5 kg	17 - 19 cm	27 - 30 cm	300000	-	300040
0	1,5 - 2,5 kg	19 - 21 cm	29 - 34 cm	300001	300021	300041
1	2,5 - 3,5 kg	21 - 24 cm	33 - 40 cm	300002	300022	300042
2	3,5 - 4,5 kg	24 - 27 cm	36 - 43 cm	300003	300023	300043
3	4,5 - 6 kg	27 - 30 cm	40 - 47 cm	300004	300024	300044
4	6 - 7,5 kg	30 - 33 cm	45 - 52 cm	300005	300025	300045
5	7,5 - 10 kg	33 - 36 cm	50 - 57 cm	300006	300026	300046
6	10 - 13 kg	36 - 39 cm	55 - 62 cm	300007	300027	300047
7	13 - 17 kg	39 - 42 cm	59 - 66 cm	300008	300028	300048
8	17 - 21 kg	42 - 46 cm	63 - 80 cm	300009	300029	300049
9	21 - 25 kg	46 - 50 cm	67 - 82 cm	300010	300030	300050
10	25 - 30 kg	50 - 54 cm	72 - 87 cm	300011	300031	300051
11	30 - 35 kg	54 - 59 cm	76 - 91 cm	300012	300032	300052
12	35 - 40 kg	59 - 64 cm	80 - 95 cm	300013	300033	300053
13	40 - 45 kg	64 - 69 cm	84 - 99 cm	300014	300034	300054
14	45 - 50 kg	69 - 74 cm	90 - 105 cm	300015	300035	300055
15	50 - 55 kg	74 - 79 cm	94 - 109 cm	300016	300036	300056

Features & Benefits

- **Post-surgical clothing** for dogs and cats
- Provides **a barrier to contamination and patient interference**
- Embedded copper **protects and promotes healthy skin**
- **In most cases replaces the cone**, allowing freedom of movement and natural behaviour during recovery
- **Easy-On Step-In Zipper** allows for a quick and less stressful application for all involved
- **Protects against the sun** (30 UVA and UVB Protection)
- **Ergonomic design** with stretch fit
- 65% cotton, 30% Polyester, 5% Elastane

How to measure



*The weight is the most accurate measure to find the right size; however, for breeds such as greyhounds and bulldogs the weight cannot be used.



@Luna&Milo

lunaandmilo.com

LUNA & MILO®

LOVING PET CARE