

# InSight<sup>®</sup> DB

## Veterinary Chemistry Analyser



- Low cost, All Inclusive Reagent, Support & Maintenance Package
- Results Transmitted via InSight p-LINK to Practice Management Systems
- Regular Preventative Maintenance Visits by our Nationwide Team of Field Service Technicians
- Out of Hours Technical Support
- Fixed Price for 5 Years



**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.

# InSight<sup>®</sup> DB Chemistry Analyser

The InSight DB Veterinary Chemistry Analyser is fully automated, uses liquid reagents, and is approximately 70% cheaper to run than dry chemistry.

## FEATURES

- User friendly intuitive veterinary software
- Onboard refrigerated reagent/sample compartment
- Minimal user maintenance
- Fast throughput - example batch times:  
9 mins 1x GHP, 13 mins 2x GHP, 16 mins 3x GHP
- Continuous loading, no need to wait for batch completion
- Stat function for urgent tests
- Fully automated with auto-reflex dilution of over range results
- Flexible testing, single parameter or up to 20 user defined profiles
- Small footprint - 414 (W) x 565 (D) mm
- Online remote assistance directly to your analyser
- Long life ceramic syringes to ensure accuracy and precision
- Environmentally friendly semi-disposable reaction cuvettes with high efficiency wash system



## CUSTOMER DEFINED PROFILES AND/OR TESTS CAN BE PROGRAMMED FROM THE FOLLOWING LIST:

Albumin, Alkaline Phosphatase, ALT, AST, Amylase, Bilirubin (Total), Calcium, Cholesterol, CK, Creatinine, CRP, Glucose, GGT, Phosphate, Protein (Total) and Urea.

BHB, Bile Acids, Chloride, Copper, Fibrinogen, Fructosamine, GLDH, GSH-Px (as Selenium), Haemoglobin, IgG, LDH, Lipase, Magnesium, NEFA, SAA, T4 (Total Thyroxin), Triglycerides, Uric Acid, Urine Creatinine & Urine Protein, Zinc, Phenobarbitone, Iron, Ferritin and TIBC.

## DIMENSIONS/WEIGHT

Dimensions: 414 (W) x 565 (D) x 385 (H) mm  
Weight: 28kg  
Water Consumption: <3L/H

**Available as a low cost package including analyser, reagents and service. Suitable for any size practice, equine practices or for multi-site placements.**



**Fully supported nationwide by our specialist team of field service technicians.**

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**